# AKRON BIO

### ABSTRACT

The discovery of the type II prokaryotic CRISPR "immune system"

Akron's optimized homogenization matrix which allows isolation of has allowed the development of an RNA-guided genome editing the Tag-less nuclease of interest at a relatively high purity after tool that is simple, easy, and quick to implement. Since then, several capture chromatography (Cation Exchange Chromatography). A approaches to genome editing have emerged. The applications of second chromatography step consisting of Size Exclusion such nucleases are ever-expanding and were recently validated Chromatography (SEC) resolves impurities (up to 100 kDa) from with the US FDA approval of Casgevy, the first gene editing therapy the 160 kDa nuclease using an enhanced buffer system. A 97.6% for the treatment of sickle cell disease. The approval of Casgevy, purity, measured by RP-HPLC, was achieved after three coupled with the many gene editing programs necessitates chromatography columns. A final chromatography step, designed versatile, cost-effective enzyme production and purification solely for the reduction of endotoxin, achieved a reduction from processes. A well-known strategy to facilitate protein purification is  $\sim$  69.8 EU/mg to levels as low as 1.0 EU/mg. The efficiency of the to engineer and express the protein with a Tag (e.g., Histidine Tag). first CEX was shown to be scalable from a 5 mL column volume (CV) This allows the use of affinity-based purification to achieve high to 5,000 mL CV. This newly developed purification process using purity. Cleavage of the Tag post-purification may be necessary, as cGMP compatible resins, chemicals, and supplies is being assessed its retention can cause interference with the protein's intended for scalability at Akron to establish a cost-effective nuclease function. In this work, we describe a purification process using purification process that does not require affinity chromatography.

# BACKGROUND

The goal of this project is to develop a scalable, cGMP compliant, pure, biologically active, 160 kDa endonuclease. Assessment of Tag-less endonuclease purification process. It will utilize each unit of operation and optimizations will be performed based homogenization, ion exchange chromatography, size exclusion upon chromatograms and analytical results. chromatography, and Tangential Flow Filtration to yield a highly

# MATERIALS & METHODS

#### Materials

- Cell pellet generated from Nuclease Production Platform
- ÄKTA avant and ÄKTAexplorer, Repligen KR2i, GEA Panda Homogenizer, Waters H Class UPLC System (Xbridge Premier Protein BEH C4)
- Clarification filters, Centrifuge
- TFF Cassettes of appropriate size to avoid breakthrough of protein of interest
- HiScale columns with applicable cGMP compatible resins (CEX, SEC, AEX)
- Denovix Spectrophotometer for concentration by UV280 nm as applicable, BioTek Spectrophotometer for quantification of material using Bradford, 4-20% SDS-PAGE for visual in-process assessment, Western Blot (SDS-PAGE, primary and secondary antibody, fluorophore, transfer buffer, PVDF membrane, milk blocking buffer, wash buffer, blotting reagent buffer), ELISA kit, Biological Activity (10X Reaction buffer, agarose gel, nuclease, sgRNA, 1 kb DNA template, thermocycler, Qiagen PCR Purification Kit), Charles River Endotoxin system and cartridges for endotoxin quantification

#### Methods

- Prepacked columns and loose cation exchange resins were obtained from industry recognized vendors and screened. All resins were sanitized pre- and post-use and were stored in appropriate storage buffers.
- The resins that supported the optimal outcome were selected.
- Buffer compositions were optimized to improve respective unit operations. Chromatography methods were modified to support scalability, resolution, and improvement in yields.
- Densitometry by SDS-PAGE was used for in-process evaluation of fraction relative purities. SDS-PAGE gels were imaged on the GelDoc and fraction purity was assessed using ImageLab 6.1 software. This generated a band % for the protein of interest in each lane, which corresponds to the purity. Average purity from fractions containing the protein of interest are reported. Error bars represent the standard error of the fractions' measured purity.
- For final purity measurement Waters H Class UPLC based in-house method was used.



# Optimized Homogenization Matrix to Produce High Purity, Tag-less, 160kDa Nuclease Using Ion Exchange and Size Exclusion Chromatography

Virginia Nash, Brittany Sangiovanni, Eva Geygan, Dylan Vasconcelles, Dennis Concepcion, Carson White, and Mahajoub Bello-Roufai Akron Bio, 600 Tallevast Road, Suite 201 & 202, Sarasota, FL 34243

Variable	Initial Homogenization Matrix	Optimized Homogenization Matrix	
Salt Concentration	Low	High	
Reducing Agent	A	В	
Surfactant	No	Yes	
pH Resuspension Buffer	7.2	7.2	
<b>Resuspension Solids %</b>	15.5%	13.4%	
Homogenization Cycles	5	5	
% Reduction OD	95.3%	94.5%	





Figure 2: Chromatograms of CEX resin screening. (A) Resin A: 7.7 mL CV. (B) Resin B: 5 mL CV. (C) Resin C: 4.7 mL C n fractions were analyzed on SDS-PAGE. The protein of interest average relative purity was assessed by SDS-PAGE-based Densitometry using ImageLab 6.1. Error bars represented standard error. No peak was eluted from Resin C under the load, wash, and elution parameters evaluated.





Theoretical Yield (%) vs. Salt Concentration Medium Salt Concentration

SEC elution fractions from each condition were analyzed on SDS-PAGE. The average purity of collected fractions for the protein of interest was assessed by SDS-PAGE based densitometry using ImageLab 6.1. Error bars represent Standard Error of purity for multiple SEC cycles (A, B, C). Theoretical yield was determined by measuring peak concentration UV at 280 nm. Error bars on theoretical yield represent Standard Error of yield obtained from multiple SEC cycles (D, E, F)

Sample	(EU/mg)	Total Endotoxin in Sample (EU)	LRV	Endotoxin Clearance
Load 11.06 mL	69.8	596.65	1942	08 37%
Pool 11.06 mL	< 0.1	< 10	1.045	70.3270
Load 280.78 mL	142	4987.04	0 0 0 0 0 0	09 4 4 9/
Pool 10X 280.78 mL	2.135	~ 67.04	0.023	70.00%



Load Concentration (mg/mL)

Table B: Data indicating the reduction of endotoxin using anion exchange chromatography







Lane	Sample
1	Linear Ladder
2	Template (Undigested)
3	Control Nuclease
4	Akron Nuclease
5	Linear Ladder

Figure 9: Agarose Gel based Nuclease activity test. The template (757 bp) was incubated in the absence of Nuclease or in the presence of an equivalent concentration of the newly purified Akron and Control Nuclease. Digest mixtures were nalyzed using Agarose Gel electrophoresis.



Figure 10: (A) UPLC chromatogram indicating a purity of 97.6% for final material. (B) SDS-PAGE of final purified material in final formulation buffer both with and without glycerol.

Conc B

Figure 1: Effect of Homogenization matrix on main elution peak purity. CEX Elution peaks obtained from Clarified lysates of both homogenization matrix were analyzed on SDS-PAGE. The average relative purity of the protein of interest for each CEX eluate was assessed by SDS-PAGE-based Densitometry using ImageLab 6.1. Error bars represent standard error for the average purity of the CEX eluate.

	1	2	3	4	5	6	7	Lane	Sample
250-	-							1	Precision Plus
							-	2	Control
150-	-	-	-	-			-	3	TFF1 Load
100- 75-	-				111		-	4	TFF1 Load 0.2µn Filtered
50-	-				-		-	5	TFFI Retentate
37-	-							6	TFF1 Permeate
25- 20-	=		-	-	-		-	7	TFF1 Recovery Flush
15-									

Figure 5: SDS-PAGE of CEX elution material ultrafiltration prior to SEC.

Lane	Salt	Load Conc. (mg/mL)	% CV
1	Low	1.46	2
2	Low	1.46	2
3	Low	1.46	2
4	Low	1.84	3
5	Low	1.84	3
6	Low	1.84	3
7	Low	1.17	3
8	Low	1.17	2
9	High	1.73	1
10	Low	1.73	1
11	Medium	1.49	1
12	Medium	1.49	1
13	Medium	1.49	1

Figure 8: Chromatogram of SEC run on a 292 mL column using optimized Salt concentration. Load concentration. and Load volume conditions. The chromatogram showed reduction of the aggregate peak that eluates before the main product peak, and resolution separating the main peak from the impurity peaks.



Lane	Sample
1	Precision Plus
2	Control
3	TFF Load
4	Retentate w/ Glycerol
5	Recovery Flush w/ Glycerol
6	UF Permeate w/Glycerol
7	UFDF Permeate w/ Glycerol
8	UFDF Permeate 2 w/ Glycerol
9	Blank
10	Control
11	Retentate no Glycerol
12	Recovery Flush no Glycerol
13	UF Permeate no Glycerol
14	UFDF Permeate w/out Glycero

#### DISCUSSION

#### Homogenization

Modification of the resuspension buffer to a higher salt concentration with a surfactant and reducing agent B aided separation of the lower molecular weight bands from the protein of interest in the CEX unit operation (Table A and Figure 1).

# CEX

Three CEX resins were screened. Resin B was selected based upon resolution, capture of protein of interest, and attainability (off shelf/minimal lead times). Modifications to residence time during elution enhanced resolution for the peak of interest (Figure 3.A.). Improvements on load duration were obtained by reducing the volume of the load adjustment, while maintaining comparable elution profiles. The CEX elution profile was transitioned from a linear gradient to a step gradient based upon conductivity, allowing scalability and reproducibility (Figure 3). Development was performed on 5 mL, 71.02 mL, and 298.4 mL columns. SDS-PAGE supported capture of the protein of interest from column volumes ranging from 5 mL to 5.09 L (Figure 4).

#### TFF

A 30 kDa regenerated cellulose cassette allowed concentration of the CEX elution material prior to application on the SEC column without permeate breakthrough (Figure 5).

# SEC

The size exclusion columns were loaded with concentrated CEX elution material. Optimizations were performed on load ratio and buffer composition (Figure 6 and Figure 7). Resolution of peaks in the SEC chromatograms improved (Figure 8), allowing a final product purity of up to 97.6% based upon HPLC (Figure 10.A.). Purified material was obtained from SEC columns ranging from 292 mL to 21.04 L.

# AEX

The AEX unit operation successfully removed endotoxin from ~69.8 EU/mg to levels as low as <0.1 EU/mg (Table B).

This process flow allowed purification of the endonuclease while maintaining biological activity (Figure 9). The undigested template shows the expected 757 bp band. The digested template was cut into 515 and 242 bp bands as anticipated. This indicates that endonuclease remains active throughout processing.

#### REFERENCES/ACKNOWLEDGEMENTS

The whole Akron Bio team has been instrumental in the production of this work.

Thank you to our vendors for supplying resin samples for screening and their support.